



# Universal Pipette Tips

*For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.*

## Universal Pipette Tips with Advanced Features and Premium Materials



### PT-10-R

#### 10 $\mu$ L Universal Tip


Material: PP (Polypropylene)  
Length: 32mm  
Capacity: 10 $\mu$ L  
Order QTY: 10 Racks per pack  
(96 tips in each rack)



### PT-200-R

#### 200 $\mu$ L Universal Tip

Material: PP (Polypropylene)  
Length: 51mm  
Capacity: 200 $\mu$ L  
Order QTY: 10 Racks per pack  
(96 tips in each rack)



### PT-1000-R

#### 1000 $\mu$ L Universal Tip

Material: PP (Polypropylene)  
Length: 78mm  
Capacity: 1000 $\mu$ L  
Order QTY: 10 Racks per pack  
(100 tips in each rack)

### Product Features:

- Produced in a Class 100,000 clean room
- Pharmaceutical grade Polypropylene
- DNase-free, RNase-free, pyrogen-free
- No endotoxin or PCR inhibitor
- Works with various pipettes such as Gilson, Eppendorf, Thermo (Finn)
- Hydrophobic filter avoids cross contamination
- Low retention tips reduce residue on the inner wall to ensure accurate results
- Gamma Irradiated Sterile

