



MINI CENTRIFUGES

These personal centrifuges are ideal for quick spins, micro-filtration of samples, cell separations, and other routine laboratory procedures.

FUN NEW COLORS!



Purple Mini Centrifuge
MLX-304A



White Mini centrifuge
MLX-306A



Grey Mini Centrifuge
MLX-310A



PRODUCT FEATURES:

- Ergonomically designed with a small footprint
- High grade plastic cover for clear visibility
- Smooth and quiet operation
- Automatic stop when lid is opened
- Digital control delivers precise speed
- Fits a variety of tubes with included accessories

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative equipment and Cryo Storage Solutions for today's scientist.

Specifications

Model	MLX-304A	MLX-306A	MLX-310A
Product Name	Purple Mini Centrifuge	White Mini Centrifuge	Grey Mini Centrifuge
Speed	4,000 RPM	6,000 RPM	10,000 RPM
Centrifugal Force	900xg	2,000xg	5,000xg
Noise Level	≤60 dB		
Capacity	6 x 1.5mL / 6x0.5mL / 16x0.2mL		
Accessories Included in Standard Packaging	Rotor Heads: MLXF-7 (6x1.5/2mL) and MLXF-3 (16x0.2mL) Adaptors: MLXF-1 (6x0.5mL) and MLXF-2 (6x0.2mL)		
Body Color	Purple	White	Grey
Operating Environment	Ambient Temperature: 5°C~40°C Relative Humidity: ≤ 80%		
Power Supply	AC110V, 60Hz		
Max. Power	8W	10W	30W
Dimensions WxDxH	155 x 140 x 115 mm (6.1 x 5.5 x 4.5 in)		
Net Weight	≤ 1kg (2.2 lb)		
Safety Feature	Stops when Lid is Opened		
Protection Level	IP21		

Accessories



MLXF-1



MLXF-2



MLXF-3



MLXF-7

Name	Part Number	PCS/Pack
Adaptors for 0.5mL Tubes	MLXF-1	6
Adaptors for 0.2mL Tubes	MLXF-2	6
Rotor Head for 16 x 0.2mL Tubes	MLXF-3	1
Rotor Head for 6 x 1.5/2mL Tubes	MLXF-7	1