




# DISPOSABLE SEROLOGICAL PIPETTES

Sterile, high quality medical grade polystyrene Disposable Serological Pipettes designed for accurate dispensing.



## PRODUCT FEATURES:

- Made of Highly Transparent medical-grade Polystyrene
- Each Pipette is Individually Wrapped, Sterile with high-quality filters
- Manufactured in Class 100,000 Cleanroom
- Includes Both Ascending and Descending Graduations that are  $\pm 2\%$  Accuracy
- Integrity Test Performed on Each Pipette
- Free of RNase, DNase, DNA, PCR inhibitors, and Pyrogens
- Color Coded by Volume

## Specifications

Part Number	SP-5-B	SP-10-B	SP-25-B	SP-50-B
Volume	5mL	10mL	25mL	50mL
Outer Dimensions	8 x 35mm	10 x 340mm	15 x 310mm	18 x 355mm
Material	Polystyrene	Polystyrene	Polystyrene	Polystyrene
Color	Blue	Orange	Purple	Black
Quantity	50/Pk	50/Pk	50/Pk	50/Pk

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative Benchtop Equipment and Cryo Storage Solutions for today's scientist.