



FOR STANDARD 2" BOXES
RACKS UF-222 TO UF-572



UF-222

SIDE ACCESS UPRIGHT FREEZER RACKS

- Designed to hold standard 2" cardboard (white or yellow), plastic (polypropylene or polycarbonate), or metal (stainless or aluminum) boxes.
- Constructed of corrosion-resistant stainless steel with double welds and smooth edges.
- Security lock devices are available for each rack.

Type of Boxes

1. WHITE 2" CARDBOARD BOX

2. POLYCARBONATE PLASTIC BOX

3. METAL BOX S-STAINLESS STEEL/A-ALUMINUM

SB2CC: 5 1/4 x 5 1/4 x 2"
133 x 133 x 51mm



1.



2.

PB2-100: 5 1/4 x 5 1/4 x 2"
133 x 133 x 52MM




3.

SB2S/SB2SXX: 5 1/4 x 5 1/4 x 2"
133 x 133 x 51MM

FOR STANDARD 2" BOXES RACKS UF-222 TO UF-572

SPECIFICATIONS

RACKS UF-222 TO UF-572



Part Number	Description	Dimension D x H x W	Compartment Size d x h x w	Capacity (Box)	Lock Device Option**	Shipping Weight lb(Kg)
UF-222	Rack Only	11 1/2 x 4 1/2 x 5 1/2 inches (292x114x140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	2 x 2 = 4	DP-1005	2.3(1.04)
UF-232	Rack Only	11 1/2 x 6 5/8 x 5 1/2 inches (292 x 169 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	2 x 3 = 6	DP-1007	3.1(1.40)
UF-242	Rack Only	11 1/2 x 8 7/8 x 5 1/2 inches (292 x 225 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	2 x 4 = 8	DP-1009	3.9(1.76)
UF-252	Rack Only	11 1/2 x 11 x 5 1/2 inches (292 x 280 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	2 x 5 = 10	DP-1011	4.5(2.06)
UF-262****	Rack Only	11 1/2 x 12 5/8 x 5 1/2 inches (292 x 321 x 140mm)	5 3/8 x 2 1/16 x 5 7/16 inches (136 x 52.6 x 138mm)	2 x 6 = 12	DP-1013	5.2(2.35)
UF-272	Rack Only	11 1/2 x 15 3/8 x 5 1/2 inches (292 x 391 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	2 x 7 = 14	DP-1015	6.0(2.72)
UF-332	Rack Only	16 7/8 x 6 6/58 x 5 1/2 inches (429 x 168 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138mm)	3 x 3 = 9	DP-1607	4.3(1.97)
UF-342 UF-342-W-XX*	Rack Only Rack w/2" box and divider *	16 7/8 x 8 7/8 x 5 1/2 inches (429 x 225 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	3 x 4 = 12	DP-1609	5.5(2.49) 5.5(2.49)+ 12w***
UF-352	Rack Only	6 7/8 x 11 x 5 1/2 inches (429 x 280 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	3 x 5 = 15	DP-1611	6.6(2.97)
UF-362****	Rack Only	16 7/8 x 12 5/8 x 5 1/2 inches (429 x 321 x 140mm)	5 3/8 x 2 1/16 x 5 7/16 inches (136 x 52.6 x 138mm)	3 x 6 = 18	DP-1613	7.5(3.38)
UF-372	Rack Only	16 7/8 x 15 3/8 x 5 1/2 inches (429 x 391 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	3 x 7 = 21	DP-1615	8.5(3.87)
UF-432	Rack Only	22 3/16 x 6 6/58 x 5 1/2 inches (429 x 168 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	4 x 3 = 12	DP-2107	5.8(2.62)
UF-442 UF-442-W-XX*	Rack Only Rack w/2" box and divider *	22 3/16 x 8 7/8 x 5 1/2 inches (429 x 225 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	4 x 4 = 16	DP-2109	7.1(3.23) 7.1 (3.23)+ 16w***
UF-452 UF-452-W-XX*	Rack Only Rack w/2" box and divider *	22 3/16 x 11 x 5 1/2 inches (564 x 280 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	4 x 5 = 20	DP-2111	8.5(3.86) 8.5(3.86)+20w***
UF-462**** UF-462-W-XX*	Rack Only Rack w/2" box and divider *	22 3/16 x 12 5/8 x 5 1/2 inches (429 x 321 x 140mm)	5 3/8 x 2 1/16 x 5 7/16 inches (136 x 52.6 x 138mm)	4 x 6 = 24	DP-2113	9.7(4.38) 9.7(4.38)+24w***
UF-472	Rack Only	22 3/16 x 15 3/8 x 5 1/2 inches (429 x 391 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	4 x 7 = 28	DP-2115	11.2(5.06)
UF-532	Rack Only	27 5/16 x 6 5/8 x 5 1/2 inches (693 x 168 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	5 x 3 = 15	DP-2707	7.0(3.18)
UF-542 UF-542-W-XX*	Rack Only Rack w/2" box and divider *	27 5/16 x 8 7/8 x 5 1/2 inches (693 x 225 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	5 x 4 = 20	DP-2709	8.6(3.91) 8.6(3.91)+20w***
UF-552	Rack Only	27 5/16 x 11 x 5 1/2 inches (693 x 280 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	5 x 5 = 25	DP-2711	10.3(4.66)
UF-562****	Rack Only	27 5/16 x 12 5/8 x 5 1/2 inches (693 x 321 x 140mm)	5 3/8 x 2 1/16 x 5 7/16 inches (136 x 52.6 x 138mm)	5 x 6 = 30	DP-2715	11.7(5.32)
UF-572	Rack Only	27 5/16 x 15 3/8 x 5 1/2 inches (693 x 391 x 140mm)	5 3/8 x 2 3/16 x 5 7/16 inches (136 x 55 x 138 mm)	5 x 7 = 35	DP-2717	13.5(6.14)

*To order, replace XX with either 81 or 100 cell divider

** Lock device is ordered separately

*** w=0.2 2lb(0.1kg) for cardboard box, 0.37lb(0.17kg) for plastic box, 0.84lb(0.38kg) for stainless steel box, 0.53lb(0.24kg) for aluminum box.

****UF-262, UF-362, and UF-562 can be fitted with 2" cardboard boxes only; plastic boxes do not fit in these racks.

