




# PIPETTE TIPS

With Advanced Feature and Premium Materials



Left to Right: PT-10-R, PT-200-R, PT-1000-R

## PRODUCT FEATURES:

- Produces in Class 100,000 Cleanroom
- Pharmaceutical Grade Polypropylene
- DNase-free, RNase-free, Pyrogen-free
- No endotoxins or PCR inhibitors
- Works with various pipettes including Gilson, Eppendorf, Thermo (Finn)
- Gamma Irradiated Sterile
- Hydrophobic filter avoids cross contamination
- Low retention tips reduce residue on the inner wall to ensure accurate results

## Specifications

	PT-10-R	PT-200-R	PT-1000-R
Material	Polypropylene	Polypropylene	Polypropylene
Length	32mm	51mm	78mm
Capacity	10µl	200µl	1000µl
Order Qty	10 racks per pack (96 tips in each rack)	10 racks per pack (96 tips in each rack)	10 racks per pack (100 tips in each rack)

*For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative Benchtop Equipment and Cryo Storage Solutions for today's scientist.*