

Rockers

◆ Description

The 2D Platform Rocker is excellent for a gentle mix of solutions, reagents, blotting and staining applications. It has a high-quality motor for durability and ultra-quiet operation. Space saving design and simple to operate make this a popular item for any laboratory.

◆ Features

- ☐ All models come with a dimpled non-skid mat.
- ☐ Large platform and hooks for rubber strips accommodate many applications.
- ☐ Units provide high capacity with rack accessory.
- ☐ All units are compatible with cold room and incubator environment.
- ☐ NS-05A , NS-14UA models have speed adjustment, continuous and timer mode with automatic stop.
- ☐ NS-14UA model provide additional feature of tilt angle adjustment from 1-15°.



2D Rocker **NS-04A**



Digital Rocker **NS-05UA**



Adjustable Rocker **NS-14UA**

◆ Specifications

Name	2D Rocker	Digital Rocker	Adjustable Rocker
Description	Fixed Speed 2D Rocker	2D Rocker w / adjustable Speed and Time	High profile 2D rocker w / adjustable time, speed and angle
Model	NS-04A	NS-05UA	NS-14UA
For regions other than North America, check the model number in the " Order Information Table "			
Speed	24rpm	2~80rpm	1~50rpm
Speed Accuracy	NA	± 1rpm	
Tilt Angle	10°		1~15° (adjustable in 1° increment)
Timer	NA	0~99h59min	1s~9999min
Platform Dimensions	307x267mm (12.1"x10.5")		
Max. Load	2kg (4.4 lb)		4.5kg (9.9 lb)
Display	NA	LCD	LED
Accessories included in Standard Packaging	NSA02, 6 Pcs of NSA01		
Ambient Temperature	5°C~40°C		
Relative Humidity	≤80%		
Power Supply	AC110V / 60Hz	AC100~240V, 50 / 60Hz	
Max. Power	4W	10W	13.5W
Dimensions (LxWxH)	355x240x123mm (14x9.5x4.8")		415x295x209mm (16.3x11.6x8.2")
Net Weight	4.5kg (9.9 lb)	6kg (13.2 lb)	14kg (30.9 lb)
Protection Level	IP21		
Certifications	NA		cCSAus, CE, RoHS, WEEE

◆ Accrssories

Part Number	Name	Specifications
NSA01	Rubber Strips	WxD: 180x7mm (7.09x0.28")
NSA02	Dimpled Rubber Mat	WxD: 298x258mm (11.73x10.16")
NSA04	Stacking Tray Assembly	WxDxH: 295x322x109mm (11.61x12.68x4.29")



◆ Order Information Table

Order Model	Power Supply	Plug Type
NS-04A	AC110V±11V, 60Hz	US
NS-05UA	AC100~240V, 50Hz / 60Hz	
NS-04E*	AC220V±22V, 50Hz	Eur
NS-05UE	AC100~240V, 50Hz / 60Hz	
NS-04EN*	AC220V±22V, 50Hz	UK
NS-05UEN	AC100~240V, 50Hz / 60Hz	
NS-04AU*	AC220V±22V, 50Hz	AU
NS-05UAU	AC100~240V, 50Hz / 60Hz	
NS-04SA*	AC220V±22V, 50Hz	SA
NS-05USA	AC100~240V, 50Hz / 60Hz	
NS-14UW	AC100~240V, 50Hz / 60Hz	Eur, UK, AU

* Note: For NS-04E, NS-04EN, NS-04AU, NS-04SA rockers, rocking speed is 20rpm.