



Centrifuge Tube Racks

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

Easy access and Storage of 15ml and 50ml Centrifuge Tubes

UFD-LT15-2
15ml (2) drawer



UFD-LT50-2
50ml (2) drawer



UFD-LT15/50
15ml (1) drawer
50ml (1) drawer



Product Features:


- Racks are specially designed to provide an ideal storage solution for large quantities of samples in 15ml and 50ml centrifuge tubes
- Drawer Racks facilitate both access to and storage of samples without having to remove the whole rack from the freezer.
- Security lock devices are available for each rack
- Made of corrosion resistant Pharmaceutical Grade Stainless Steel
- Available in 15ml, 50ml and 15/50ml drawers

Ordering Information:

Content to Store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension L x H x W	Capacity	Lock Device Option *	Shipping Weight (lb/kg)
	16 1/2 inch 419 mm	4 15/16 inch 126 mm	1	UFD-T15-1	Rack with Drawer for 15ml Centrifuge Tubes	16 11/16 x 4 9/16 x 5 1/2 inch 424 x 126 x 140 mm	60 Tubes (15 ml)	DD-04	4.0 / 1.83
	22 inch 559 mm			UFD-LT15-1		22 3/16 x 4 15/16 x 5 1/2 inch 564 x 126 x 140 mm	80 Tubes (15 ml)		5.2 / 2.35
	26 3/4 inch 679 mm			UFD-XLT15-1		26 15/16 x 4 9/16 x 5 1/2 inch 685 x 126 x 140 mm	104 Tubes (15 ml)		6.1 / 2.75
	16 1/2 inch 419 mm	9 15/16 inch 252 mm	2	UFD-T15-2		16 11/16 x 9 9/16 x 5 1/2 inch 424 x 252 x 140 mm	120 Tubes (15 ml)	DD-09	7.4 / 3.34
	22 inch 559 mm			UFD-LT15-2		22 3/16 x 9 15/16 x 5 1/2 inch 564 x 252 x 140 mm	160 Tubes (15 ml)		9.3 / 4.22
	26 3/4 inch 679 mm			UFD-XLT15-2		26 15/16 x 9 15/16 x 5 1/2 inch 684 x 252 x 140 mm	208 Tubes (15 ml)		10.8 / 4.90
	16 1/2 inch 419 mm	4 15/16 inch 126 mm	1	UFD-T50-1	Rack with Drawer for 50ml Centrifuge Tubes	16 11/16 x 4 9/16 x 5 1/2 inch 424 x 126 x 140 mm	30 Tubes (50 ml)	DD-04	3.9 / 1.78
	22 inch 559 mm			UFD-LT50-1		22 3/16 x 4 15/16 x 5 1/2 inch 564 x 126 x 140 mm	39 Tubes (50 ml)		5.1 / 2.30
	26 3/4 inch 679 mm			UFD-XLT50-1		26 15/16 x 4 9/16 x 5 1/2 inch 685 x 126 x 140 mm	51 Tubes (50 ml)		5.9 / 2.68
	16 1/2 inch 419 mm	9 15/16 inch 252 mm	2	UFD-T50-2		16 11/16 x 9 9/16 x 5 1/2 inch 424 x 252 x 140 mm	60 Tubes (50 ml)	DD-09	7.2 / 3.25
	22 inch 559 mm			UFD-LT50-2		22 3/16 x 9 15/16 x 5 1/2 inch 564 x 252 x 140 mm	78 Tubes (50 ml)		9.1 / 4.12
	26 3/4 inch 679 mm			UFD-XLT50-2		26 15/16 x 9 15/16 x 5 1/2 inch 684 x 252 x 140 mm	102 Tubes (50 ml)		10.4 / 4.76
	16 1/2 inch 419 mm	4 15/16 inch 126 mm	2	UFD-T15/50	Rack with Drawer for 15/50ml Centrifuge Tubes	16 11/16 x 9 15/16 x 5 1/2 inch 424 x 252 x 140mm	60 Tubes (15 ml) 30 Tubes (50 ml)	DD-04	4.6 / 2.10
	22 inch 559 mm			UFD-LT15/50		22 3/16 x 9 15/16 x 5 1/2 inch 564 x 252 x 140mm	80 Tubes (15 ml) 39 Tubes (50 ml)		5.5 / 2.50
	26 3/4 inch 679 mm			UFD-XLT15/50		26 15/16 x 9 15/16 x 5 1/2 inch 685 x 252 x 140mm	104 Tubes (15 ml) 51 Tubes (50 ml)		5.5 / 2.50



(1) Drawer model



Dimensions