



# Cryovials

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

## Crystal Clear Leak Proof Cryovials for Biological Sample Storage at Temperatures Down to -196°C

**IN STOCK NOW**



### Product Features:

- Caps and Vials are molded from Medical Grade Polypropylene ensuring the optimum seal to prevent leakage
- The Internal Thread Screw Cap design & O-ring reduces contamination
- Sterile (Gamma Radiation) and free of endotoxin, RNase, DNase, DNA, PCR inhibitors, and pyrogens
- Rounded internal bottom of vial promotes complete pour out
- Crystal Clear Transparent vial for easy observation
- Large white writing area allows easy marking
- Manufactured in Cleanroom environment
- Designed for freezing samples in any application

Note: Not to be in the liquid phase of liquid nitrogen tanks

### Specifications:

| Part Number | Material | Outer Dimension | Volume Capacity | Temperature Range | Qty. in PK | Description                    |
|-------------|----------|-----------------|-----------------|-------------------|------------|--------------------------------|
| CVI-2.0-B   | PP       | Ø 12.5 x 50mm   | 2.0ml           | -196 ~ 121°C      | 50         | Internal Thread, Self Standing |