




# Multi-Position Magnetic Stirrers

For over 25 years Crystal Technology & Industries has been driven to deliver exceptional value to the scientific marketplace. Our engineers collaborate with top-tier research institutions to design & develop innovative Benchtop Equipment and Cryo Storage Solutions for today's Scientist.

This Multi-Position Magnetic Stirrer base is fully submersible, Ideal for Water Bath or C02 Incubator. Remote Panel allows for full control while in operation



## Product Features:

- Sealed stainless steel plate protects against intrusion and guarantees reliability.
- Waterproof body for use in water bath. Protection Level: IP21 (Control), IP68S (Platform)
- Magnetic coil drive mechanism provides long & quiet operation.
- Five stirring levels provide excellent mixing capability, including liquids with high viscosity.
- Previous setting is automatically saved when restarting the unit.
- See [crystalindustries.com](http://crystalindustries.com) for full listing of accessories available to fit your application.

## Product Specifications:

Product	Multi-Position Magnetic Stirrers		
Model	MS-04BSUA	MS-06BSUA	MS-15BSUA
Speed	80 ~ 2000 rpm		
Speed Accuracy	± 10%		
Max. Length for Stir Bar	50mm (1.97")		20mm (.79")
Platform Dimensions	250 x 254mm (9.84 x 10")	384 x 254 (15.12 x 10")	385 x 254mm (15.16 x 10")
Capacity	1000ml		100ml
Display	LCD		
Viscosity Levels	5		
Number of Stirring Positions	4	6	15
Accessories Included	1x (One) Stirrer Bar MSA04		1x (One) Stirrer Bar MSA06
Dimensions (WxDxH)	260 x 264 x 43mm (10.24 x 10.39 x 1.69")		390 x 264 x 42mm (15.35 x 10.39 x 1.65")
Net Weight	5.5kg (12.13 lb)		9kg (19.84 lb)