




CRYOVIALS

Crystal Clear Leak Proof Cryovial for Biological Sample Storage at Temperatures Down to 196°C



PRODUCT FEATURES:

- Internally threaded cap
- Manufactured in Cleanroom environment (Class 100,000)
- Free of RNase, DNase, DNA, PCR inhibitors, and Pyrogens
- Sterile (Gamma Irradiated)
- Caps and Vials molded from Medical Grade Polypropylene ensuring optimum seal to prevent leakage
- Designed for freezing samples in various application including vapor phase Liquid Nitrogen

(Not designed to be used in the liquid phase of liquid nitrogen)

IN STOCK NOW

Specifications

Part Number	Material	Outer Dimesions	Volume Capacity	Temperature Range	Qty in PK	Description
CVI-2.0-B	Polypropylene	12.5 x 50mm	2.0ml	-196 ~ 121°C	Pk 50	Internal Thread, Self Standing

For over 25 years Crystal Technology and Industries has been driven to deliver exceptional value to the scientific marketplace by developing innovative Benchtop Equipment and Cryo Storage Solutions for today's scientist.