





**CRYSTAL TECHNOLOGY
& INDUSTRIES, INC.**

**Freezer Racks & Boxes
for Ultra-Low Freezers,
Freezers, and Refrigerators**



2025 Cryo Storage Solutions

Our Most Popular Premium Racks & Boxes



Table of Contents

Standard Side Access Racks for 2 & 3" Boxes	4
Adjustable Freezer Racks	5
Sliding Tray Freezer Racks for 2 & 3" Boxes	6
Premium Drawer Freezer Racks for 2 & 3" Boxes	7
Vertical Racks for 2 & 3" Boxes	8
Specialty Vertical Racks for 2" Boxes	9
Vertical Racks for SBS Plates, Racks/Boxes.	10
Premium Drawer Freezer Racks for 15mL & 50mL Centrifuge Tubes	11
Side Access Racks for 15mL & 50mL Centrifuge Tube Boxes	11
Standard Side Access Racks for SBS Plates, 96 Well & 384 Well Plates.	12
Premium Drawer Racks for 96-Well & 384 Well Plates	12
Freezer Boxes	13 & 14
Rack Configurator	15

Crystal offers the widest range of Freezer Racks and Boxes to fit most brands/models of Ultra-Low Freezers, Upright Freezers, Chest Freezers and Liquid Nitrogen Tanks. This Mini Catalog highlights our most popular models, please see our website crystalindustries.com or ask your representative for more sizes and configurations.

All Crystal racks:

- Are constructed of high grade, corrosion resistant stainless steel
- Come with an unconditional lifetime warranty on workmanship
- Manufactured using the latest precision technology, robotic welding, CNC machines, and CAD design software
- Feature double welds on all seams
- Standard racks for 2", 3", 4" and 5" Boxes, Microplates, Specialty Tube Racks and 15mL & 50mL Centrifuge Tubes
- Most Racks in stock for quick delivery
- Cryo Storage Solutions for LN2 Tanks available, contact us!



Don't see what you need? Contact us for easy, economical Custom Racks!

We make choosing the right Freezer Racks and Boxes to optimize your storage space easy, check out our website for our Freezer Rack Configurator. crystalindustries.com/freezer-rack-configurator/

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Standard Side Access Racks for 2" & 3" Boxes



◆ Features

- Constructed of high grade corrosion resistant stainless steel.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Compartments accept all cardboard, plastic and metal boxes.
- Handles located on both ends for easy handling.

For 2" Freezer Boxes

Item#	Description / Capacity	L x H x W
UF-332	3 x 3 Rack / for 9 Boxes	16-7/8 x 6-5/8 x 5-1/2"
UF-342	3 x 4 Rack / for 12 Boxes	16-7/8 x 8-7/8 x 5-1/2"
UF-352	3 x 5 Rack / for 15 Boxes	16-7/8 x 11 x 5-1/2"
UF-432	4 x 3 Rack / for 12 Boxes	22-3/16 x 6-5/8 x 5-1/2"
UF-442	4 x 4 Rack / for 16 Boxes	22-3/16 x 8-7/8 x 5-1/2"
UF-452	4 x 5 Rack / for 20 Boxes	22-3/16 x 11 x 5-1/2"
UF-462*	4 x 6 Rack / for 24 Boxes	22-3/16 x 12-5/8 x 5-1/2"
UF-472	4 x 7 Rack / for 28 Boxes	22-3/16 x 15-3/8 x 5-1/2"
UF-542	5 x 4 Rack / for 20 Boxes	27-5/16 x 8-7/8 x 5-1/2 "
UF-552	5 x 5 Rack / for 25 Boxes	27-5/16 x 11 x 5-1/2"

For 3" Freezer Boxes

Item#	Description / Capacity	L x H x W
UF-323	3 x 2 Rack / for 6 Boxes	16-7/8 x 6-5/8 x 5-1/2"
UF-333	3 x 3 Rack / for 9 Boxes	16-7/8 x 9-3/8 x 5-1/2"
UF-423	4 x 2 Rack / for 8 Boxes	22-3/16 x 6-5/8 x 5-1/2"
UF-433	4 x 3 Rack / for 12 Boxes	22-3/16 x 9-3/18 x 5-1/2"
UF-443	4 x 4 Rack / for 16 Boxes	22-3/16 x 12-5/8 x 5-1/2"
UF-453	4 x 5 Rack / for 20 Boxes	22-3/16 x 15-9/16 x 5-1/2"
UF-523	5 x 2 Rack / for 10 Boxes	27-5/16 x 6-5/8 x 5-1/2"
UF-533	5 x 3 Rack / for 15 Boxes	27-5/16 x 9-1/2 x 5-1/2"

*UF-462 is designed to hold standard 2" cardboard boxes only, plastic boxes do not fit.

Ask your representative about security lock accessories available for these racks. We have over 1,000+ sizes & styles of racks. Ask your representative if you don't see the size you need.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Adjustable Freezer Racks for 2" & 3" Boxes



◆ Features

- Allows you to store a variety of different boxes in one rack.
- Adjustable & removable stainless steel clips to accommodate a wide variety of box heights.
- Lighter... So it's easier to take in and out of your freezer.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFMC-33-2-3	3 x 3 Adj. Rack / for 9 Boxes	16-15/16 x 6-5/8 x 5-1/2"
UFMC-34-2-3	3 x 4 Adj. Rack / for 12 Boxes	16-15/16 x 8-7/8 x 5-1/2"
UFMC-35-2-3	3 x 5 Adj. Rack / for 15 Boxes	16-15/16 x 11 x 5-1/2"
UFMC-43-2-3	4 x 3 Adj. Rack / for 12 Boxes	22-3/8 x 6-5/8 x 5-1/2"
UFMC-44-2-3	4 x 4 Adj. Rack / for 16 Boxes	22-3/8 x 8-7/8 x 5-1/2"
UFMC-45-2-3	4 x 5 Adj. Rack / for 20 Boxes	22-3/8 x 11 x 5-1/2"
UFMC-46-2-3	4 x 6 Adj. Rack / for 24 Boxes	22-3/8 x 12-5/8 x 5-1/2"
UFMC-54-2-3	5 x 4 Adj. Rack / for 20 Boxes	27-5/16 x 8-7/8 x 5-1/2"
UFMC-55-2-3	5 x 5 Adj. Rack / for 25 Boxes	27-5/16 x 11 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFMC-33-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 6-5/8 x 5-1/2"
UFMC-34-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 8-7/8 x 5-1/2"
UFMC-35-2-3	3 x 2 Adj. Rack / for 6 Boxes	16-15/16 x 11 x 5-1/2"
UFMC-43-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 6-5/8 x 5-1/2"
UFMC-44-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 8-7/8 x 5-1/2"
UFMC-45-2-3	4 x 2 Adj. Rack / for 8 Boxes	22-3/8 x 11 x 5-1/2"
UFMC-46-2-3	4 x 2 Adj. Rack / for 12 Boxes	22-3/8 x 12-5/8 x 5-1/2"
UFMC-54-2-3	5 x 2 Adj. Rack / for 10 Boxes	27-5/16 x 8-7/8 x 5-1/2"
UFMC-55-2-3	5 x 2 Adj. Rack / for 10 Boxes	27-5/16 x 11 x 5-1/2"

*All Above Adjustable Racks Come Complete with Clips.

Additional Clips		
Item#	Description / Capacity	L x H x W
UFMC-L	"L" clip	3-5/8 x 3/8 x 1"
UFMC-T	"T" clip	3-5/8 x 3/8 x 1/2"

Need help with which rack fits your freezer?

Not storing your samples in boxes?
Ask your representative about our new bin style racks.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Sliding Tray Freezer Racks for 2" & 3" Boxes



◆ Features

- Convenience of accessing each tray individually.
- Premium performance at a value price.
- Space saving design to maximize your storage needs.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFS-332	3 x 3 Slide Rack / for 9 Boxes	16-7/8 x 6-5/8 x 5-1/2"
UFS-342	3 x 4 Slide Rack / for 12 Boxes	16-7/8 x 8-7/8 x 5-1/2"
UFS-352	3 x 5 Slide Rack / for 15 Boxes	16-7/8 x 11 x 5-1/2"
UFS-432	4 x 3 Slide Rack / for 12 Boxes	22-3/8 x 6-5/8 x 5-1/2"
UFS-442	4 x 4 Slide Rack / for 16 Boxes	22-3/8 x 8-7/8 x 5-1/2"
UFS-452	4 x 5 Slide Rack / for 20 Boxes	22-3/8 x 11 x 5-1/2"
UFS-522	5 x 2 Slide Rack / for 10 Boxes	27-1/16 x 4-9/16 x 5-1/2"
UFS-532	5 x 3 Slide Rack / for 15 Boxes	27-1/16 x 6-5/8 x 5-1/2"
UFS-542	5 x 4 Slide Rack / for 20 Boxes	27-1/16 x 8-7/8 x 5-1/2"
UFS-552	5 x 5 Slide Rack / for 25 Boxes	27-1/16 x 11 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFS-333	3 x 3 Slide Rack / for 9 Boxes	16-7/8 x 9-1/2 x 5-1/2"
UFS-343	3 x 4 Slide Rack / for 12 Boxes	16-7/8 x 12-11/16 x 5 1/2"
UFS-353	3 x 5 Slide Rack / for 15 Boxes	16-7/8 x 15-13/16 x 5 1/2"
UFS-433	4 x 3 Slide Rack / for 12 Boxes	22-3/8 x 9-1/2 x 5-1/2"
UFS-443	4 x 4 Slide Rack / for 16 Boxes	22-3/8 x 12-11/16 x 5-1/2"
UFS-453	4 x 5 Slide Rack / for 20 Boxes	22 3/8 x 15-13/16 x 5 1/2"
UFS-533	5 x 3 Slide Rack / for 15 Boxes	27-1/16 x 9-1/2 x 5-1/2"
UFS-543	5 x 4 Slide Rack / for 20 Boxes	22-1/16 x 12-11/16 x 5-1/2"

Cardboard boxes are recommended for use with the sliding Tray Racks. Looking for a smaller or larger rack? Ask your representative, as we have additional sizes.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Premium Drawer Freezer Racks for 2" & 3" Boxes



◆ Features

- Fully framed providing even greater sample protection.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Compartment options to accept all cardboard, plastic and metal boxes.
- Label place card holders for easy identification.

For 2" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFD-332	3 x 3 Drawer Rack / for 9 Boxes	16-3/4 x 7-1/8 x 5-1/2"
UFD-342	3 x 4 Drawer Rack / for 12 Boxes	16-3/4 x 9-7/16 x 5-1/2"
UFD-352	3 x 5 Drawer Rack / for 15 Boxes	16-3/4 x 11-13/16 x 5-1/2"
UFD-432	4 x 3 Drawer Rack / for 12 Boxes	22-5/16 x 7-1/8 x 5-1/2"
UFD-442	4 x 4 Drawer Rack / for 16 Boxes	22-5/16 x 9-7/16 x 5-1/2"
UFD-452	4 x 5 Drawer Rack / for 20 Boxes	22-5/16 x 11-13/16 x 5-1/2"
UFD-462	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 14-1/8 x 5-1/2"
UFD-462-S*	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 12-5/8 x 5-1/2"
UFD-462-SA**	4 x 6 Drawer Rack / for 24 Boxes	22-5/16 x 13-1/4 x 5-1/2"
UFD-542	5 x 4 Drawer Rack / for 20 Boxes	27 x 9-7/16 x 5-1/2"
UFD-552	5 x 5 Drawer Rack / for 25 Boxes	27 x 11-13/16 x 5-1/2"

For 3" Freezer Boxes		
Item#	Description / Capacity	L x H x W
UFD-323	3 x 2 Drawer Rack / for 6 Boxes	16-3/4 x 6-7/8 x 5-1/2"
UFD-333	3 x 3 Drawer Rack / for 9 Boxes	16-3/4 x 10-1/16 x 5-1/2"
UFD-423	4 x 2 Drawer Rack / for 8 Boxes	22-5/16 x 6-7/16 x 5-1/2"
UFD-433	4 x 3 Drawer Rack / for 12 Boxes	22-5/16 x 10-1/16 x 5-1/2"
UFD-443	4 x 4 Drawer Rack / for 16 Boxes	22-5/16 x 13-1/2 x 5-1/2"
UFD-453	4 x 5 Drawer Rack / for 20 Boxes	22-5/16 x 16-9/16 x 5-1/2"
UFD-523	5 x 2 Drawer Rack / for 20 Boxes	27-1/16 x 6-7/8 x 5-1/2"
UFD-533	5 x 3 Drawer Rack / for 15 Boxes	27-1/16 x 10-1/16 x 5-1/2"

Constructed of high grade corrosion resistant stainless steel with precision craftsmanship. Ask your representative about label cards and security lock accessories for these racks.

*Top-Loading for Panasonic Freezers for 2" Cardboard Boxes only

**Side-Loading for Panasonic Freezers for 2" Polycarbonate Boxes

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Vertical Racks for 2" & 3" Boxes



◆ Features

- Security lock devices are available for each rack.
- Constructed of corrosion resistant stainless steel.
- Locking Rods included for each rack.

For 2" Freezer Boxes, Stainless Steel Rack

Item#	Description / Capacity	L x H x W
CF-4-2	4 Boxes	9-1/8 x 5-5/8 x 5-1/2"
CF-5-2	5 Boxes	11-5/16 x 5-5/8 x 5-1/2"
CF-6-2	6 Boxes	13-1/2 x 5-5/8 x 5-1/2"
CF-7-2	7 Boxes	15-11/16 x 5-5/8 x 5-1/2"
CF-8-2	8 Boxes	17-7/8 x 5-5/8 x 5-1/2"
CF-9-2	9 Boxes	20-1/16 x 5-5/8 x 5-1/2"
CF-10-2	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
CF-11-2	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
CF-12-2	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
CF-13-2	13 Boxes	28-1/2 x 5-5/8 x 5-1/2"
CF-14-2	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
CF-15-2	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

For 3" Freezer Boxes, Stainless Steel Rack

Item#	Description / Capacity	L x H x W
CF-4-3	4 Boxes	13 x 5-5/8 x 5-1/2"
CF-5-3	5 Boxes	16-3/16 x 5-5/8 x 5-1/2"
CF-6-3	6 Boxes	19-3/8 x 5-5/8 x 5-1/2"
CF-7-3	7 Boxes	22-1/2 x 5-5/8 x 5-1/2"
CF-8-3	8 Boxes	25-11/16 x 5-5/8 x 5-1/2"
CF-9-3	9 Boxes	28-1/2 x 5-5/8 x 5-1/2"
CF-10-3	10 Boxes	32-1/16 x 5-5/8 x 5-1/2"

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Specialty Vertical Racks for 2" Boxes

◆ Features

- Made of aluminum.
- Security lock devices are available for each rack.



For 2" Freezer Boxes, Aluminum Rack

Item#	Description / Capacity	L x H x W
CFA-10-2	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
CFA-11-2	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
CFA-12-2	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
CFA-13-2	13 Boxes	28-9/16 x 5-5/8 x 5-1/2"
CFA-14-2	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
CFA-15-2	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

*Other sizes available

◆ Features

- Made of corrosion resistant stainless steel.
- Security lock devices are available for each rack.
- Spring Clips securely hold boxes in rack compartments without using locking rods and allow for easy removal of boxes.



For 2" Freezer Boxes, Stainless Steel Rack with Spring Clip

Item#	Description / Capacity	L x H x W
CF-9-2-S	9 Boxes	20-1/16 x 5-5/8 x 5-1/2"
CF-10-2-S	10 Boxes	22-5/16 x 5-5/8 x 5-1/2"
CF-11-2-S	11 Boxes	24-1/2 x 5-5/8 x 5-1/2"
CF-12-2-S	12 Boxes	26-5/8 x 5-5/8 x 5-1/2"
CF-13-2-S	13 Boxes	28-1/2 x 5-5/8 x 5-1/2"
CF-14-2-S	14 Boxes	31-1/16 x 5-5/8 x 5-1/2"
CF-15-2-S	15 Boxes	33-1/4 x 5-5/8 x 5-1/2"

◆ Features

- Made of corrosion resistant stainless steel.
- Racks are ideal for light load and easy handing of sample storage in freezers or liquid nitrogen tanks.



For 2" Freezer Boxes, Vertical Half Racks

Item#	Description / Capacity	L x H x W
LR-3-2	3 Boxes	16-1/2 x 2-1/4 x 5-9/16"
LR-4-2	4 Boxes	21-5/16 x 2-1/4 x 5-9/16"
LR-5-2	5 Boxes	26-3/4 x 2-1/4 x 5-9/16"
LR-6-2	6 Boxes	32-13/16 x 2-1/4 x 5-9/16"

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Vertical Racks for SBS Plates, Racks/Boxes



For 15mL & 50mL
Centrifuge Tube Boxes



For 96-Deep
Well Plates



For Specialty
Tube Racks



For 100 Cell Hinged
Plastic Boxes

CFLB Series For 15mL & 50mL Centrifuge Boxes

Item#	Description / Capacity	L x H x W
CFLB-4	4 Boxes	20-3/8 x 6 x 6-3/16"
CFLB-5	5 Boxes	25-3/8 x 6 x 6-3/16"
CFLB-6	6 Boxes	30-3/8 x 6 x 6-3/16"

Other sizes and configurations available see our website or ask your representative.

CFMX Series For Specialty Tube Racks

Item#	Description / Capacity	L x H x W
CFMX-9-2	9 Boxes	22-5/8 x 3-5/8 x 5-1/2"
CFMX-10-2	10 Boxes	25-1/8 x 3-5/8 x 5-1/2"
CFMX-11-2	11 Boxes	27-5/8 x 3-5/8 x 5-1/2"
CFMX-12-2	12 Boxes	30-1/16 x 3-5/8 x 5-1/2"
CFMX-13-2	13 Boxes	32-9/16 x 3-5/8 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

CFDP Series For 96-Deep Well Plates

Item#	Description / Capacity	L x H x W
CFDP-9	9 Boxes	16-5/8 x 3-5/8 x 5-1/2"
CFDP-10	10 Boxes	18-7/16 x 3-5/8 x 5-1/2"
CFDP-11	11 Boxes	20-1/4 x 3-5/8 x 5-1/2"
CFDP-12	12 Boxes	22-1/8 x 3-5/8 x 5-1/2"
CFDP-13	13 Boxes	23-15/16 x 3-5/8 x 5-1/2"
CFDP-14	14 Boxes	25-3/4 x 3-5/8 x 5-1/2"
CFDP-15	15 Boxes	27-9/16 x 3-5/8 x 5-1/2"
CFDP-16	16 Boxes	29-7/16 x 3-5/8 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

CFHT Series For 100 Cell Hinged Plastic Boxes

Item#	Description / Capacity	L x H x W
CFHT-8	8 Boxes	18-3/4 x 5-7/8 x 6-1/4"
CFHT-9	9 Boxes	21 x 5-7/8 x 6-1/4"
CFHT-10	10 Boxes	23-3/8 x 5-7/8 x 6-1/4"
CFHT-11	11 Boxes	25-5/8 x 5-7/8 x 6-1/4"
CFHT-12	12 Boxes	27-15/16 x 5-7/8 x 6-1/4"
CFHT-13	13 Boxes	30-1/4 x 5-7/8 x 6-1/4"

Other sizes and configurations available see our website or ask your representative.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Premium Drawer Freezer Racks for 15mL & 50mL Centrifuge Tubes

◆ Features

- Available in a wide variety of formats:
1 or 2 drawer, 15mL, 50mL & 15mL & 50mL.
- Precision craftsmanship with tolerances of +/- 1 mm.
- Easy drawer access to each tube and no need for boxes.



For 15mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
UFD-T15-2	120 Tubes (15mL)	16-11/16 x 9-15/16 x 5-1/2"
UFD-LT15-2	160 Tubes (15mL)	22-3/16 x 9-15/16 x 5-1/2"
UFD-XLT15-2	200 Tubes (15mL)	26-15/16 x 9-15/16 x 5-1/2"

For 50mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
UFD-T50-2	60 Tubes (50mL)	16-11/16 x 9-15/16 x 5-1/2"
UFD-LT50-2	78 Tubes (50mL)	22-3/16 x 9-15/16 x 5-1/2"
UFD-XLT50-2	102 Tubes (50mL)	26-15/16 x 9-15/16 x 5-1/2"

For 15mL & 50mL Centrifuge Tubes

Item#	Description / Capacity	L x H x W
UFD-T15/50	60 Tubes (15mL) + 30 Tubes (50mL)	16-11/16 x 9-15/16 x 5-1/2"
UFD-LT15/50	80 Tubes (15mL) + 39 Tubes (50mL)	22-3/16 x 9-15/16 x 5-1/2"
UFD-XLT15/50	100 Tubes (15mL) + 51 Tubes (50mL)	26-15/16 x 9-15/16 x 5-1/2"

Storing a different tube size? We have more racks for those too... Ask your representative.

Side Access Racks for 15mL & 50mL Centrifuge Tube Boxes



For 15mL & 50mL Centrifuge Tube Boxes

Item#	Description / Capacity	L x H x W
UFLB-22	2 x 2 Rack / for 4 Boxes	12-3/4 x 10-3/16 x 5-1/2"
UFLB-32	3 x 2 Rack / for 6 Boxes	18-13/16 x 10-3/16 x 5-1/2"
UFLB-42	4 x 2 Rack / for 8 Boxes	24-13/16 x 10-3/16 x 5-1/2"

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Standard Side Access Racks for SBS Plates, 96 Well & 384 Well Plates

◆ Features

- Locking rods are included.
- Security lock devices are available for each rack.
- Made of corrosion resistant stainless steel.



Item#	Description / Capacity	L x H x W
UFMP-509L	70 Plates w / lids 80 Plates w / o lids	18-13/16 x 9-3/16 x 5-1/2"
UFMP-511L	85 Plates w / lids 95 Plates w / o lids	18-13/16 x 11 x 5-1/2"
UFMP-609L	84 Plates w / lids 96 Plates w / o lids	22-1/2 x 9-3/16 x 5-1/2"
UFMP-611L	96 Plates w / lids 120 Plates w / o lids	22-1/2 x 11 x 5-1/2"
UFMP-709L	98 Plates w / lids 112 Plates w / o lids	26-1/8 x 9-3/16 x 5-1/2"
UFMP-711L	112 Plates w / lids 140 Plates w / o lids	26-1/8 x 11 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

Premium Drawer Racks for 96-Well & 384-Well Plates

◆ Features

- Fully framed providing greater sample protection.
- Compartment fits all standard 96-Well & 384-Well Plates
- Label place card holders for easy identification




Item#	Description / Capacity	L x H x W
UFDMP-54	4 drawers 60 Plates w / lids 80 Plates w / o lids	18-3/16 x 9-3/4 x 5-1/2"
UFDMP-55	5 drawers 75 Plates w / lids 100 Plates w / o lids	18-3/16 x 12-1/16 x 5-1/2"
UFDMP-64	4 drawers 72 Plates w / lids 96 Plates w / o lids	21-3/4 x 9-3/4 x 5-1/2"
UFDMP-65	5 drawers 90 Plates w / lids 120 Plates w / o lids	21-3/4 x 12-1/8 x 5-1/2"
UFDMP-66	6 drawers 108 Plates w / lids 144 Plates w / o lids	21-3/4 x 14-1/2 x 5-1/2"
UFDMP-74	4 drawers 84 Plates w / lids 112 Plates w / o lids	25-1/4 x 9-3/4 x 5-1/2"

Other sizes and configurations available see our website or ask your representative.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Freezer Boxes



Water Proof Cardboard Boxes

Item#	Description / Capacity	L x W x H	Color
SB2CC-H-81	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
SB2CC-H-81-P-W	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	Purple
SB2CC-H-81-G-W	2" Water Proof Cardboard w/ Hinged Lid w/81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	Green

Water Repellent Cardboard Boxes

Item#	Description / Capacity	L x W x H	Color
SB2CC	Standard 2" Cardboard Box	5-1/4 x 5-1/4 x 2"	White
SB2CC-81	Standard 2" Cardboard Box, with 81 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
SB2CC-100	Standard 2" Cardboard Box, with 100 Place Cell Divider inserted	5-1/4 x 5-1/4 x 2"	White
SB2CC-LP*	Low Profile Water Repellent 2" Cardboard Box	5-5/16 x 5-5/16 x 1-15/16"	White
SB2C*	Standard 2" Cardboard Box	5-1/4 x 5-1/4 x 2"	Yellow
SB3C*	Standard 3" Cardboard Box	5-1/4 x 5-1/4 x 3"	White
SB3.75C*	Standard 4" Cardboard Box	5-1/4 x 5-1/4 x 3-7/8"	White
SB4C*	Standard 5" Cardboard Box	5-13/16 x 5-13/16 x 4-15/16"	White

*Available with cell dividers


Standard Box Divider

Item#	Description / Capacity	L x D x W	Cell Opening
CD-64	8 x 8 = 64 Cells	4-7/8 x 4-7/8 x 1-5/16"	9/16 inch
CD-81-W	9 x 9 = 81 Cells	4-7/8 x 4-7/8 x 1-5/16"	1/2 inch
CD-100	10 x 10 = 100 Cells	4-7/8 x 4-7/8 x 1-5/16"	7/16 inch
CD-144	12 x 12 = 144 Cells	4-7/8 x 4-7/8 x 1-5/16"	3/8 inch

Many more sizes available see our website or contact your representative

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Freezer Boxes



◆ Features

Stainless steel boxes come with or without cardboard dividers and fit into standard freezer racks.

- 2" and 3" boxes available
- With or without cell dividers
- Fits into standard freezer boxes
- Drain holes



◆ Features

Aluminum boxes come with or without cardboard dividers and fit into standard freezer racks.

- Standard 2" Aluminum box with drain holes
- Standard 3" Aluminum box with drain holes
- Cell dividers available for each
- Fits into standard freezer boxes

PolyCarbonate, Polypropylene and Plastic Boxes

Item#	Description / Capacity	L x H x W
PB2-100*	10 x 10 Polycarbonate Box	5-1/4 x 5-1/4 x 2-1/16"
PB2-81*	9 x 9 Polycarbonate Box	5-1/4 x 5-1/4 x 2-1/16"
PB4-81*	9 x 9 Polycarbonate Box	5-1/4 x 5-1/4 x 3-11/16"
MPB2-25*	5 x 5 Polycarbonate Mini Box	3 x 3 x 2-1/16"
PPP2-100	10 x 10 Polypropylene Box	5-1/4 x 5-1/4 x 2"
PPP2-81	9 x 9 Polypropylene Box	5-1/4 x 5-1/4 x 2"
PP2-100	100-Cell Hinged Top Plastic Box pk1 (Natural)	5-9/16 x 6 x 2-1/8"
PP2-100-A	100-Cell Hinged Top Box pk 5 (5 colors)	5-9/16 x 6 x 2-1/8"
PP2-100X*	100-Cell Hinged Top Box (Specify 5 color: B-Blue, G-Green, O-Orange, P-Pink, Y-Yellow)	5-9/16 x 6 x 2-1/8"

*Suitable for use to -196°C

Stainless Steel Boxes

Item#	Description / Capacity	Dimension Wt×Wt×Ht (Lid)	Dimension Wb×Wb×Hb
SB2S	Standard 2" Stainless Steel Box w/Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB2S-XX*	Standard 2" Stainless Steel Box w/Divider & Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB3S	Standard 3" Stainless Steel Box w/Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm
SB3S-XX*	Standard 3" Stainless Steel Box w/Divider & Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm

* To order, replace XX with 144, 100, 81, 64, 49, 36, 25 or 16 for cell dividers.

Aluminum Boxes

Item#	Description / Capacity	Dimension Wt×Wt×H	Dimension Wb×Wb×Hb
SB2A	Standard 2" Aluminum Box w/Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB2A-XX*	Standard 2" Aluminum Box w/Divider & Drain Hole	5-1/4×5-1/4×2" 133×133×51mm	5 x 5 x 1-15/16" 127.7×127.7×49mm
SB3A	Standard 3" Aluminum Box w/Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm
SB3A-XX*	Standard 3" Aluminum Box w/Divider & Drain Hole	5-1/4×5-1/4×3" 133×133×75mm	5 x 5 x 2-15/16" 127.7×127.7×74mm


* To order, replace XX with 144, 100, 81, 64, 49, 36, 25 or 16 for cell dividers.

Freezer Racks & Boxes for Ultra-Low Freezers, Freezers, and Refrigerators

Rack Configurator

We make choosing the right Freezer Racks and Boxes to optimize your storage space easy, check out our website for our Freezer Rack Configurator.
crystalindustries.com/freezer-rack-configurator/

The screenshot shows the Crystal Industries website navigation bar with the logo on the left, a search bar, and links for Products, Distributors, Freezer Rack & Box Configurator (circled in red), and Contact. Below the navigation bar is a blue header with the text "Freezer Rack & Box Configurator". The main content area is divided into two columns. The left column is titled "Freezer Rack Configurator" and includes a note: "Note: If your freezer Make and/or Model is not shown, choose 'Not on List' and you will be prompted to enter your freezer dimensions. The configurator will then automatically determine which freezer rack will fit your freezer appropriately." Below the note are five dropdown menus: "Freezer Type", "Freezer Make", "Freezer Model", "Box Type", and "Rack Type". The right column is titled "Freezer Box Configurator" and includes a note: "Note: If your Tube Type or Tube Size is not shown, choose 'Not on List' and you will be prompted to enter your Tube dimensions. The configurator will then automatically determine which Boxes will fit your Tubes appropriately." Below the note are three dropdown menus: "Tube Type", "Tube Size", and "Box Type".



Crystal Industries Sales & Product Support:

1600 S. Airport Dr. | McKinney, TX 75069

972-934-2525

crystalindustries.com

