




Centrifuge Tube Freezer Racks for Ultra-Low Freezers, Freezers, Refrigerators, and Benchtops

Easy access and storage of 15mL and 50mL Centrifuge Tubes



- Racks are specially designed to provide an ideal storage solution for large quantities of samples of 15mL and 50mL centrifuge tubes
- Drawer Racks facilitate both access to and storage of samples without having to remove the whole rack from the freezer
- Sizes available to fit most brands of Ultra-Low Freezers, Freezers and Refrigerators.
- Security lock devices are available for each rack
- Made of corrosion resistant Pharmaceutical Grade Stainless Steel
- Available in 15mL, 50mL and 15/50mL drawers



UFD-XLT15-2
2 Drawer
Freezer Rack
15mL Tubes

UFD-T50-2
2 Drawer
Freezer Rack
50mL Tubes

UFD-LT15/50
2 Drawer Freezer Rack
15mL/50mL Tubes

Content to store	Rack Length	Rack Height	Number of Drawers	Part Number	Description	Dimension LxHxW	Capacity	Lock Device Option*
15mL Tubes	16-1/2 in 419 mm	4-15/16 in 126 mm	1	UFD-T15-1	Rack with drawer for 15mL Centrifuge Tubes	16-11/16 x 4-15/16 x 5-1/2 in 424 x 126 x 140 mm	60 tubes	DD-04
	22 in 559 mm			UFD-LT15-1		22-3/16 x 4-15/16 x 5-1/2 in 564 x 126 x 140 mm	80 tubes	
	26-3/4 in 679 mm			UFD-XLT15-1		26-15/16 x 4-15/16 x 5-1/2 in 685 x 126 x 140 mm	100 tubes	
	26-3/4 in 679 mm	9-15/16 in 252 mm	2	UFD-XLT15-2		26-15/16 x 9-15/16 x 5-1/2 in 684 x 252 x 140 mm	200 tubes	DD-09
50mL Tubes	16-1/2 in 419 mm	4-15/16 in 126 mm	1	UFD-T50-1	Rack with drawer for 50mL Centrifuge Tubes	16-11/16 x 4-15/16 x 5-1/2 in 424 x 126 x 140 mm	30 tubes	DD-04
	22 in 559 mm			UFD-LT50-1		22-3/16 x 4-15/16 x 5-1/2 in 564 x 126 x 140 mm	39 tubes	
	26-3/4 in 679 mm			UFD-XLT50-1		26-15/16 x 4-15/16 x 5-1/2 in 685 x 126 x 140 mm	51 tubes	
	16-1/2 in 419 mm	9-15/16 in 252 mm	2	UFD-T50-2		16-11/16 x 9-15/16 x 5-1/2 in 424 x 252 x 140 mm	60 tubes	DD-09
	26-3/4 in 679 mm			UFD-XLT50.2		26-15/16 x 9-15/16 x 5-1/2 in 684 x 252 x 140 mm	102 tubes	
15/50mL Tubes	16-1/2 in 419 mm	4-15/16 in 126 mm	2	UFD-T15/50	Rack with drawer for 15/50mL Centrifuge Tubes	16-11/16 x 9-15/16 x 5-1/2 in 424 x 252 x 140 mm	60 tubes (15mL) 30 tubes (50mL)	DD-04
	22 in 559 mm			UFD-LT15/50		22-3/16 x 9-15/16 x 5-1/2 in 564 x 252 x 140 mm	80 tubes (15mL) 39 tubes (50mL)	
	9-15/16 in 252mm			UFD-XLT15/50		26-15/16 x 9-15/16 x 5-1/2 in 684 x 252 x 140 mm	100 tubes (15mL) 51 tubes (50mL)	

Contact your Crystal Industries Sales Representative for more information on these and our full line of Cryo Storage Solutions.

