

IS-18CA Large Capacity & IS-ACMA Compact CO2/NonCO2 and 170L Space saving CO2/NonCO2 Incubator Shaker in Stock with 2 Year Warranty!



- Models that Feature CO2, Shaking, Humidity and Refrigeration
- Designed to Provide the Optimum Environment for Samples
- Reliable, Quiet and Built to Last
- Large Capacity 6 x 5L & Compact Models 4 x 5L with Stackable Design up to 3 Units
- Space Saving 170L CO2/NonCO2 Models with Static Shelf for Parallel Cell Cultivation and can be stacked 2 units high
- Models feature Touchscreen Technology, Pull-out Platform, Built-in UV Sterilization, Microplate Holder Platform
- Choose from 3mm, 25mm or 50mm Orbit Models
- Wide Range of Accessories Available including Clamps, Sticky Mats and Tube Holders

Top Hinge Incubator Shaker



- Models Available with Cooling
- 4 x 2L Capacity
- 25mm Orbit, 30-400 RPM
- Clear Tempered Glass for Easy Viewing
- Wide Range of Accessories Available

IS-RSDAA Bench-Top Incubator Shaker



- Versatile platform with a wide selection of accessories
- Precise temperature control (ambient +5°C~ 80°C)
- Speeds of up to 300 RPM with an orbit diameter of 26mm (1.02")
- The Sliding platform feature provides easy sample handling with Max Load 5x1L
- Fits conveniently under countertops

CI-191CXA CO2 Incubator



- Intuitive User Interface with 7" Touchscreen
- Precise IR CO2 Sensor
- Easy High Heat Sterilization
- Temperature Controlled Door for Higher Uniformity
- HEPA Filter provides Class 100 Purification
- 3 Year Warranty

NEW! CO2 Resistant Orbital Shaker for Existing CO2 Incubators



- Designed for use in CO2, Standard or Refrigerated Incubators
- Magnetic Orbital Drive System for precise speed control, (no belts!) 19mm Orbit, 30-300 RPM
- Remote Control attaches externally to CO2 Incubator; LED display allows users to set and see actual parameters without opening door

Crystal Industries Sales & Product Support:

Linda Hurley / lhurley@crystalindustries.com / 484-886-9655

Essential Benchtop Equipment in Stock with 2 Year Replacement Warranty!



Revolvers



Rotators



Mini-Vortexers



2D Rockers



Dry Baths



Vortex Mixers



3D Rockers



Mini-Rollers



Orbital Shakers



Mini-Centrifuges



Mag Stirrers



Multiposition Stirrers



Heating/Cooling Dry Baths



- Designed and Manufactured for Reliability and Durability with Advanced Features
- Variety of Accessories Available to Accommodate Most Sample Holders
- Models in Stock
- 2 Year Warranty

FREEZER RACKS with Lifetime Warranty & Rack Configurator to Make Selection Easy! Economical and Quick Custom Rack Process!



Side Access Racks



Drawer Style Racks



Blood Bag Cannisters



Vertical Racks



Centrifuge Tube Racks



Blood Bag Cannister Frames

- Widest Range of Racks to fit most Models of Upright & Chest Freezers and Liquid Nitrogen Tanks
- Manufactured using the latest technology, Robotic Welding, CNC Machines and CAD Design Software from High Quality Stainless Steel
- Double Welds on all Seams
- Blood Bag, Stem Cell Storage Systems
- Easy, Economical Custom Rack Process
- Most Racks in Stock for Quick Delivery

Crystal Industries Sales & Product Support:

Linda Hurley / lhurley@crystalindustries.com / 484-886-9655

